

# Efficient Capture of Suspended Tumor Cells with a Novel Immunocytochemical Microfluidic Device

Mary Dickson<sup>†, ‡</sup>, Pavel Tsinberg<sup>†</sup>, Zhongliang Tang<sup>†</sup>, Timothy Wilson<sup>†</sup>, Edward Leonard<sup>†</sup>

<sup>†</sup>Biocept, Inc., San Diego, CA

<sup>‡</sup>Department of Chemical Engineering, Columbia University, New York, NY

## Introduction

Capture of rare cells such as:


1. Circulating tumor cells (CTCs) in the blood of cancer patients
2. Fetal cells in the blood of expectant mothers

and their subsequent genetic investigation is crucial for the developing of individualized diagnoses & treatments.

A highly efficient, high purity separation of these cells (as rare as 1:10<sup>9</sup>) has been achieved by a method developed in the authors' lab. Here investigated are the optimal flow parameters for capture.


## Methodology

1. Tag CTCs with precise ratio of biotinylated antibodies




2. Flow blood and CTCs through streptavidin coated micro-channel lined with ~9,000 posts.

- Biotin – streptavidin bond only adheres CTCs to post



3. CTCs are located via fluorescent microscopy & cytogenetically interrogated in-situ (FISH, etc.)



## Results

Capture efficiency, C: The number of CTCs captured vs. the total number of CTCs.

$$C = \frac{n_{cap}}{n_0}$$

C depends on surface antigen density (x), & can be characterized as:

$$y = \frac{x}{a + x}$$




Figure 4: Series: Trials with and without other blood (bystander) cells. C is greater than 70% for antigen densities above 30,000. C is about 50% at antigen densities of 10,000.

How does antigen density affect cell binding?

- 1) Chance of adhesion

- 2) Strength of cell-post bond

\*Both scale with number of biotin ligands per cell

Probability of Capture, k: Likelihood of a cell to capture within a given row. In a given row (r) the change in number of cells in suspension (n) is proportional to n.

$$k = \frac{\ln[1-C]}{r}$$





Figure 6: Parameter: Distance along channel, each segment is 50 rows.

Left: High antigen density, k values decrease along channel at low flow rates.



Right: Low antigen density, k values are relatively constant

How does flow rate affect cell binding?

- 1) Higher flow-rates:

- Increased shear rate  
Displacement of cells from posts?
- Increased speed in higher flow-rate regimes  
Not enough time for binding to occur?

- 2) Low flow-rates:


- Streamlines develop which are farther from posts  
Some cells will never make contact?

## Conclusions

Extremely rare cells can be separated with high efficiency and high purity for in-situ inspection. This will prove useful in

- Advanced diagnostics
- Individualized treatments
- Targeted cell removal

Higher antigen densities yield higher capture. Middle flow rates yield constant k values, avoiding streamline effect illustrated below.



## References

- Liotta, L. A., Kleinerman, J., Sidel, G. M. "Quantitative Relationships of Tumor Cells, Tumor Vessels, and Pulmonary Metastases following Tumor Implantation." *Cancer Research* 34 (1974): 997-1004.
- Racila, E., Euhus, D., Weiss, A. J., Rao, C., McConnell, J., Terstappen, L. W. M. M., Uhr, J. W. "Detection and Characterization of carcinoma cells in the blood." *Proceedings of the National Academy of Sciences, USA*, 95 (1998): 4589-4594.
- Meng, S., Tripathy, D., Frenkel, E. P., Shete, S., Natfalis, E. Z., Huah, J. F., Beisch, P. D., Leitch, M., Hoover, S., Euhus, D., Hakey, B., Morrison, L., Fleming, T. P., Heryn, D., Terstappen, L. W. M. M., Fehn, T., Tucker, T. F., Lane, N., Wang, J. "Circulating Tumor Cells in Patients with Breast Cancer Dormancy." *Clinical Cancer Research* 10 (2004): 8152-8162.
- Pantel, K., Beakenhoff, R. H., Brandt, B. "Detection, clinical relevance, and specific biological properties of disseminating tumour cells." *Nature Reviews Cancer* 8 (2008): 329-340.
- Rothdorff, S., Fritsche, H., Muller, V., Rau, T., Schindlbeck, C., Rack, B., Janin, W., Coih, C., Beck, K., Janitsch, F., Jackson, S., Gomet, I., Cristofanilli, M., Pantel, K. "Detection of Circulating Tumor Cells in Peripheral Blood of Patients." *Clinical Cancer Research* 13.3 (2007): 920-928.
- Sesqui, L. V., Nagrath, S., Toner, M., Haber, D., Lynch, T. J. "The CTC-Chip: An Exciting New Tool to Detect Circulating Tumor Cells in Lung Cancer Patients." *Journal of Thoracic Oncology*, 4.3 (2009): 281-283.
- Handy, B. C. "Conference Scene: Advancing toward personalized medicine." *Biomarkers Medicine* 4.1 (2010): 129-131.
- Tsinberg, P., Tang, Z. Device for cell separation and analysis and method of using. United States Patent 7,695,956 B2, 2011.

## Acknowledgements

Many thanks to Lina Ebio for her skilled technical work, as well as my advisors: Tom Burns & Kathe Houghtaling

## More information?

Biocept, Inc., San Diego, CA.  
[www.biocept.com](http://www.biocept.com)

Mary Dickson  
Columbia University Department  
of Chemical Engineering  
[mary.dickson@gmail.com](mailto:mary.dickson@gmail.com)